

RECOMMENDED VENDORS FOR CONSUMABLES AND REAGENTS

EU Version

Consumables

NAME	DESCRIPTION	BRAND	PRODUCT NUMBER
Plastic Crimp/Snap vials	Crimp/Snap polypropylene 11 mm vials. Volume: 800 µL	ThermoFisher Scientific	C4011-11
Glass vials	1.5 ml Snap ring vials. Dimensions: 32 x 11.6 mm (clear), wide opening	BGB Analytik	11SR
Snap cap lids	11 mm Snap caps with 1 mm thickness. PTFE/Silicon/PTFE-septum. Red colored	Infochroma (email for quote)	G003-SCR-FKSKFK10
Nunc MicroWell 96-well plates	500 µl 96 well, V-bottom plate, polypropylene, clear. Working volume: 300µL	Sigma	P7116-1CS
96-Well deep plates	STARLAB 96. Deep-well plate, Square wells, volume = 2.2 mL. U-shaped bottom, clear	Life System Design	4Ti-0136
	Jet Biofil. Deep-well Multiwell Plates, round 1.0 mL, 96 wells, non-sterile. 24 box /96 cases, U bottom	Axonlab	RWP102096
96-Well adhesive foils	96-well Microplate foil. 100 pieces of self-adhesive, transparent plastic foils for polystyrene and polypropylene microplates	Cytiva	28975816

Reagents

NAME	DESCRIPTION	BRAND	PRODUCT NUMBER
Amine coupling kit, standard	Contains buffers and chemicals for 40 – 60 activations: <ul style="list-style-type: none"> • EDC, solid, 1 g • MES/NHS buffer solution, 50 ml • Ethanolamine* HCl, 50 ml • Borate buffer 50 ml (not indicated in the website) 	Xantec (email for order)	K AN-50
Borate buffer	Borate buffer solution 50 mL and 500 mL	Xantec (email for order, not on website)	B BELU-50 B BELU-500

